



Retina Lens

90D 78D 20D

Deeper Observation, Better Diagnosis

Specifications

90D

78D

20D



Magnification ratio	0.76 x	0.93 x	3.13 x
Field Of View(°)	74.2° ~ 89°	81.4° ~ 97°	46°
PrismWorking Distance (mm)	7.0	8.0	50.0
Clear Aperture(mm)	Ø21	Ø26	Ø48
Dimension(mm)	Ø27 x 13	Ø33 x 15	Ø56 x 20

NMPA

Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423 Fax: +86-21-54260425

Email: marketing@mediworks.biz
international@mediworks.biz

Follow us



www.mediworks.biz