




## Portable Slit Lamp S150

All in one, All in hand

# Specification

|                       |                                   |
|-----------------------|-----------------------------------|
| <b>Model No.</b>      | <b>S150</b>                       |
| Magnification         | 6 x                               |
| Filter                | Cobalt blue filter                |
| Brightness adjustment | 4 levels                          |
| Slit Length           | ≥ 12 mm                           |
| Slit Width            | 0 ~ 12 mm continuously adjustable |
| Power                 | Li-ion battery 3.7 V / 3400 mAh   |
| Illumination          | LED bulb 3 V / 1 W                |
| Working Time          | More than 6 hours                 |
| Recharging Time       | 3.5 hours                         |
| Net Weigh             | 240 g (including battery)         |

Image Gallery:



Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423 Fax: +86-21-54260425

Email: [marketing@mediworks.biz](mailto:marketing@mediworks.biz)  
[international@mediworks.biz](mailto:international@mediworks.biz)

Follow us



[www.mediworks.biz](http://www.mediworks.biz)

202306 \*Subject to change in design or specification without advance notification

**Details Make the Difference**