




**Applanation Tonometer**

**T170**

Gold standard, accurate measurement





## Accurate Measurement

With years of tests and clinical applications, Goldmann has been proved to be the gold standard for tonometry measurement. Based on the proven evidence, MediWorks T170 tonometer applied Goldmann principle in the design and manufacture.

With precision manufacturing process, MediWorks applanation tonometer offers accurate and consistent measurements which give a clear picture of patients' IOP situation to the doctors.



## Excellent Mechanics

Full attentions are paid to every detail in the manufacture of MediWorks tonometers. Workers assemble every component carefully and check tight tolerances in weight to make sure the tonometer deliver consistent measurements.



## Premium Optics

The prism on the tonometer is the key to accurate measurement. In MediWorks, our skillful workers make these prisms with great care and follow strict quality standards. We offer 3 pcs of prisms in standard package making MediWorks tonometers durable in long term use.

## Accessories



Tonometer Adaptor (Z type)  
(For Z type applanation tonometer)

(Available for S260s/ S260)

Z type



# Specifications

T170	
Range of Measurement	0 mmHg ~ 80 mmHg (0 kPa ~ 10.64 kPa)
Moving Range of the Rings	0 ~ 3 mm
Prism Diameter	7 mm
Moving Range of Prism	± 1.5 mm
Dimension	171 mm x 80.5 mm x 143.5 mm
Net Weight	284 g



Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423

Email: [marketing@mediworks.biz](mailto:marketing@mediworks.biz)  
[international@mediworks.biz](mailto:international@mediworks.biz)

Follow us



[www.mediworks.biz](http://www.mediworks.biz)