

Optometry Lift Table

ES-18A

Specification



Table top size	550 mm x 450 mm
Table elevator height	660 mm ~ 860 mm
Input	100 - 220 Vac + 10% (50/60Hz)
Fuse	5A
Weight	21KG
Maximum load	50KG



ES-40L

Specification



Table top size	1100 mm x 575 mm
Table elevator height	665 mm ~ 865 mm
Input	DC 100~220 V $\pm$ 10%, (50 Hz / 60 Hz)
Fuse	220 VA
Weight	35 KG
Maximum load	180 KG

Shanghai MediWorks Precision Instruments Co.,Ltd.

Add: Building 7, Ming Pu Plaza, No. 3279, San Lu Rd, Min Hang District, Shanghai, 201100, China

Tel: +86-21-54260421 54260423

Email: marketing@mediworks.biz
international@mediworks.biz

Follow us



www.mediworks.biz